

Lab Equipment Activity

Name:

Date:


Block:

Station 1: Using A Scale



1. What do you think a scale is used to measure? _____
2. In Canada, the SI unit **grams (g)** is used to measure mass. In order to change the unit on your scale to the one you want, you will need to press the '**M**' or the '**Unit**' button on the scale.
3. List all of the unit measurements that the scale can read.

4. Ensure that the scale displays '**g**' as your unit.
5. **Before** you measure the mass of the provided object, what *should* your scale read?
6. Hit the '**tare/zero**' button on your scale so that it reads **0.00 g**.
7. Pick up the provided object and place it on the scale. The mass of the object is _____ .
8. At times, certain chemicals in the form of a powder or a liquid require the use of a container in order to measure its mass.
9. To measure the mass of a powder, you must use a weigh boat to hold the powdered chemical. Place the weigh boat on the scale and hit the '**tare/zero**' button so that it reads **0.00 g**. Place two scoops of the chemical into the weigh boat. The mass of the chemical is _____ .
10. Remove the weigh boat and pour the chemical back into its original container.
11. To measure the mass of a liquid, we will use an empty beaker to hold the liquid. Place an empty beaker onto the scale and press the '**tare/zero**' button so that the scale reads **0.00 g**.
12. Using a **graduated cylinder**, measure out 80.0 mL of the liquid. Pour the liquid into the empty beaker. The mass of the liquid is _____ .
13. Remove the beaker and pour the liquid back into its original container.
14. Turn off the scale.



Station 2: Taking Qualitative Observations

*****Safety:** For this station, be sure that safety goggles are on AT ALL TIMES. Be sure to avoid contact with chemical solutions with eyes and all body tissues.

1. Put on a pair of **safety goggles**.
2. For each of the liquids in the beakers, make **2 qualitative observations**:

Liquid A	Liquid B	Liquid C
•	•	•
•	•	•

3. Using the **eye dropper**, place 3 – 5 drops of each liquid into separate test tubes. (Make sure that eye droppers are placed back in the original beaker)
4. Place 1 drop of Bromthymol Blue into each of the test tubes.
5. Write down 2 **qualitative observations** for each of the following reactions:

Liquid A	Liquid B	Liquid C
•	•	•
•	•	•

6. Empty the liquids in the test tubes into the sink with lots of running water. Place the test tubes upside down on a paper towel to dry.

Station 3: Taking Quantitative Observations

*****Safety:** For this station, be sure that safety goggles are on AT ALL TIMES. Be sure to avoid contact with chemical solutions with eyes and all body tissues.

1. Put on a pair of **safety goggles**.
2. Using a **graduated cylinder**, measure out EXACTLY 50.0 mL of cold water. Then pour the water into the labelled beaker.
3. Using a **graduated cylinder**, measure out EXACTLY 50.0 mL of hot water. Then pour the water into the labelled beaker.
4. Using the materials provided, make 3 **quantitative observations** for each liquid.

Cold Water	Hot Water
•	•
•	•
•	•

5. Pour the hot and cold water together into an empty beaker.
6. Make 2 more **quantitative observations** for the resultant mixture.

Resultant Mixture
•
•

7. Pour out the liquids down the sink.
8. Remove safety goggles.

Station 4: Identifying Lab Equipment

Using the provided equipment on the table, match the equipment with the names below:

- | | |
|---------------------------|-----------------------------|
| _____ 1. Beaker | _____ 8. Graduated cylinder |
| _____ 2. Thermometer | _____ 9. Eyedropper |
| _____ 3. Erlenmeyer flask | _____ 10. Safety glasses |
| _____ 4. Hot plate | _____ 11. Bunsen burner |
| _____ 5. Test tube | _____ 12. Scoopula |
| _____ 6. Test tube rack | _____ 13. Stir rod |
| _____ 7. Scale | _____ 14. Funnel |

Identify which of the following equipment is used for each of the following scenarios:

- To protect your eyes. _____
- To measure the temperature of a liquid. _____
- For approximate measurement of a liquid. _____
- For more accurate measurement of a liquid _____
- To measure the mass of a substance. _____
- To hold a test tube. _____
- To transfer small amounts of liquid from one container to another. _____
- To stir liquids. _____
- To scoop solids. _____
- To transfer liquids into a container with a small opening _____